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- [1] Babaian R, et al. NCCN Prostate Cancer Early Detection. Clinic Practice Guidelines in Oncology—V.1. National Comprehensive Cancer Network, 2004
- [2] Fleisher M, Dnistrian AM, Sturgeon CM, et al. Practice guidelines and recommendations for use tumor markers in the clinic. Washington DC: AACC Press, 2002
- [3] Semjonow A, Albrecht W, Bialk P, et al. Tumor markers in prostate cancer—EGTM recommendations. *Anticancer Res*, 1999, 19: 2799-801
- [4] Wang MC, Valenzuela IA, Murphy GP, et al. Purification of a human prostate specific antigen. *Invest Urol*, 1979, 17: 159-63
- [5] Duffy MJ, McGing P, McSweeney J. Guidelines for the Use of Tumor Markers, Produced on behalf of the Scientific Committee of the ACBI, Second edition, September 2000
- [6] Oesterling JE, Jacobsen SJ, Chute CG, et al. Serum prostate-specific antigen in a community-based population of healthy men—Establishment of age-specific reference ranges. *JAMA*, 1993, 270: 860-864
- [7] Paul R, Breul J, Hartung R. Prostate-specific antigen density and age-specific prostate-specific antigen values: the solution of prostate cancer screening *Eur Urol*, 1995, 27: 286-291
- [8] Schmid HP: Prostate specific antigen doubling time in diagnosis and follow-up of patients with prostate cancer. *Tumour Marker Update*, 1996, 8: 71-77
- [9] Amico AV, Chen MH, Roehl KA, et al. Preoperative PSA velocity and the risk of death from prostate cancer after radical prostatectomy. *N Engl J Med*, 2004, 351(2): 125-135
- [10] Bangma CH, Kranse R, Blijenberg BG. The value of screening tests in the detection of prostate cancer. Part II: Retrospective analysis of free/total prostate-specific analysis ratio, age-specific reference ranges, and PSA density. *Urology*, 1995, 46: 779-784
- [11] Catalona WJ, Smith DS, Wolfert RL, et al. Evaluation of percentage of free serum prostate-specific antigen to improve specificity of prostate cancer screening. *JAMA*, 1995, 274: 1214-1220
- [12] Catalona WJ, Partin AW, Slawin KM, et al. Use of the percentage of free prostate-specific antigen to enhance differentiation of prostate cancer from benign prostatic disease. *JAMA*, 1998, 279: 1542-1547
- [13] Collins M. Prostate Cancer: Staging of Prostate Cancer, Contemporary Issues in Prostate Cancer: A Nursing Perspective, 2002
- [14] Brawer MK. How to use PSA in the early detection or screening for prostate carcinoma. *CA Cancer J Clin*, 1995, 45: 148-64
- [15] Thomas L. Clinical laboratory diagnostics. Frankfurt Germany: TH-Books Verlagsgesellschaft mbH. 1998
- [16] Sokoll LJ, Bruzek DJ, Dua R, et al. Short-term stability of the molecular forms of prostate-specific antigen and effect on percent complexed prostate-specific antigen and percent free prostate-specific antigen. *Urology* 2002; 60: 24-30
- [17] Paus E, Nilsson O, Bormer OP. Stability of free and total prostate specific antigen in serum from patients with prostate carcinoma and benign hyperplasia. *J Urol*, 1998, 159: 1599-1605

- [18] Sturgeon C, Dati F, Duffy MJ, et al. Quality requirements and control—EGTM recommendations. *Anticancer Res* 1999, 19:2791-4
- [19] NCCLS. Primary Reference Preparations Used to Standardize Calibration of Immunochemical Assays for Serum Prostate Specific Antigen (PSA); Approved Guideline. NCCLS document I/LA19-A (ISBN 1-56238-323-X). NCCLS, Pennsylvania 19087, 1997
- [20] Catharine M. Sturgeon and Eleftherios Diamandis. Use of Tumor Markers in Testicular, Prostate, Colorectal, Breast, and Ovarian Cancers. *Laboratory Medicine Practice Guidelines*. The American Association for Clinical Chemistry, 2009
- [21] 中华医学会泌尿外科学分会. 前列腺癌诊断治疗指南, *中华现代外科学杂志*, 2006(22): 1839-1856
- [22] Carter HB, Albertsen PC, Barry MJ, et al. Early Detection of Prostate Cancer: AUA Guideline. *J Urol*. 2013 May 6. Pii: S0022-5347(13)04308-5
- [23] 那彦群、叶章群、孙颖浩, 等. *中国泌尿外科疾病诊断治疗指南(2014版)*, 北京: 人民卫生出版社, 2013.12, P.62-64
-